



The 7000C Dynamic Headspace module is a plug and play option for CDS 7000C concentrator. This module, when installed with 7000C concentrator, adds dynamic headspace capability to PAL RTC/RSI rail. The setup is capable of performing dynamic headspace analysis in a closed-system, which is USP 467, EPA 5035A and 8260C soil compliant.



User Interface of PAL DHS Plugin



**High Temp Isothermal Zones**

1. Max 300 °C Sample Vial Heater
2. Max 300 °C Valve Oven
3. Max 350 °C Transferline

**Seamless Upgrade**

1. Plug and Play Module
2. Easy Mount on PAL System
3. TCP/IP Based

**Build-in Magnetic Stirrer**

1. Select from Three Speeds
2. 8260C Compliant in Soil Mode

