

JSB GAS BAG SAMPLING MODULE

DESCRIPTION

Module for sampling from gas sample bags on the PAL3 platform. There are slots for 4 bags per module. There is no limit on the number of modules present on one PAL other than space on the x-axis of the PAL.

Suitable bags are available in various materials in various sizes from SKC:

<https://www.skcltd.com/products2/gas-sample-bags.html>

A **video** showcasing the module can be found here:

https://youtu.be/OkLtS_edpQ4


JSB PARTS:

- Gas Bag Sampling module(s) for up to 4 bags
- Needle guide adapter for syringe tool
- USB stick with passive PAL modules for import and user manual

ADDITIONAL PARTS REQUIRED:

- PAL
- Suitable syringe tool with syringe (Headspace/ITEX)


PART NUMBER	DESCRIPTION
JSBPAL 000100	Gas Bag Sampling Module bundle for up to 4 bags. Includes 1 x Gas Sampling Module, 1 x Syringe adapter, User Manual, usb stick. Let us know what software platform you are using.
JSBPAL 000101	Gas Bag Sampling Module for up to 4 bags
JSBPAL 000102	Syringe Adapter for Gas Bag Sampling Module



Parts included.



Adapter for syringe tool.



The gas sampling bag module.



One 4 bag module mounted on a PAL.

JSB United Kingdom and Ireland
 Unit 25 Graylands Estate Langhurst Wood
 Road Horsham West Sussex RH12 4QD
 t: +44 (0) 1628 854565
 e: info@go-jsb.co.uk

JSB Netherlands
 Apollweg 2B
 8239 DA Lelystad
 t: +31(0)32 087 00 18
 e: info@go-jsb.nl

www.go-jsb.com

SALES & SERVICE:

-  The Netherlands
-  Germany
-  U.K. & Ireland